Heritage Medical Products, Inc.

General Maintenance For New and Older Exam Beds

At SONOBED we strive to design our products for years of dependable service. In an effort to provide a safe product and prevent accidents or injury to your patients we suggest the following inspections be performed semi-annually, or more often under heavy patient use.

· Check the condition of the upholstery.

· Check that the plugs for the lift column, actuator(s), and hand control are fully inserted into the ports on the control box.

· Check all electrical cords for fraying or damage.

· Check the movement of all actuators and the lift column.

· Check the connections of all actuators to the frame assemblies.

· Check the connection of the middle section to the lift column.

· Check the action of the foot section. For manual foot sections, check the condition of the latch assembly and the brace.

· Check the action of the Drop Section. For manual Drop Sections, check the condition of the latch assembly and the brace.

· Check the action of the Wedge Section(s). For manual Wedge Sections, check the condition of the pivot points and the brace.

· Check the welded frame assemblies. Look for cracks or other signs of fatigue.

· Check the function of all casters. Make sure they roll smoothly and securely lock.

· Check that all bolted connections are secure.

· Check that all accessories are secured to the exam bed.